

## Product Datasheet

### H3 (phospho-Ser28) antibody (orb34390)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Rabbit polyclonal antibody to H3 (phospho-Ser28). The protein belongs to the histone H3 family, and ...   |
| <b>Species/Host</b>        | Rabbit  |
| <b>Reactivity</b>          | Human   |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Tested Applications</b> | WB  |
| <b>Immunogen</b>           | Synthetic Peptide   |
| <b>Target</b>              | H3C1  |
| <b>Form/Appearance</b>     | Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.  |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles  |
| <b>Note</b>                | For research use only   |
| <b>Isotype</b>             | Rabbit IgG  |
| <b>Clonality</b>           | Polyclonal  |
| <b>MW</b>                  | 15404   |
| <b>Uniprot ID</b>          | <a href="#">P68431</a>  |
| <b>NCBI</b>                | <a href="#">NP_003520.1</a> , <a href="#">NP_003524.1</a> , <a href="#">NP_003521.2</a> , <a href="#">NP_003526.1</a> , <a href="#">NP_066298.1</a> , <a href="#">NP_003523.1</a> , <a href="#">NP_003522.1</a> , <a href="#">NP_003528.1</a> , <a href="#">NP_003525.1</a> , <a href="#">NP_003527.1</a> |
| <b>Dilution Range</b>      | WB: 1:1000  |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |

