




Product Datasheet

HER2 (phospho-Tyr1248) antibody (orb34444)

Description	Rabbit polyclonal antibody to HER2 (phospho-Tyr1248)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT, WB
Immunogen	Synthetic Peptide
Target	ERBB2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is first purified by protein G affinity chromatography. Then, the antibody fraction is peptide affinity purified in a 2-step procedure with control and phosphorylated peptides. The phospho-specific antibody is eluted with high and low pH buffers and neutralized immediately, followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	137910
Uniprot ID	P04626
NCBI	NP_004439.2, NP_001005862.1
Dilution Range	WB: 1:8000, DOT: 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of A431 cells using...



Western blot analysis of A431 cells using...