

Product Datasheet

FGFR (phospho-Tyr766) antibody (orb34447)

Description Peptide affinity purified rabbit polyclonal antibody.

Species/Host Rabbit

Reactivity Human, Mouse

Conjugation Unconjugated

Tested Applications IHC-P, WB

Immunogen Synthetic Peptide

Target FGFR1

Form/Appearance Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal



MW 91868

Uniprot ID [P11362](#)

NCBI [NP_001167538.1](#), [NP_001167536.1](#), [NP_075593.1](#), [NP_056934.2](#), [NP_075598.2](#), [NP_001167535.1](#), [NP_001167534.1](#), [NP_075594.1](#), [NP_001167537.1](#)

Dilution Range WB: 1:1000, IHC-P: 1:50-100

Expiration Date 12 months from date of receipt.



Western blot analysis of in mouse cerebe...