
Product Datasheet

Lysozyme antibody (orb344575)

Description

Lysozyme antibody

Species/Host

Rabbit

Reactivity

Gallus

Conjugation

Unconjugated

Tested Applications

DOT, ELISA, IHC, IP, WB

Immunogen

This purified antibody was prepared from whole rabbit serum produced by repeated immunizations with full length protein corresponding to amino acids 1-129 of Hen Egg White Lysozyme.

Preservatives

0.01% (w/v) Sodium Azide

Form/Appearance

Liquid (sterile filtered)

Concentration

1.0 mg/mL

Storage

Store vial at -20° C or below prior to opening. This vial contains a relatively low volume of reagent (25 µL). To minimize loss of volume dilute 1:10 by adding 225 µL of the buffer stated above directly to the vial. Recap, mix thoroughly and briefly centrifuge to collect the volume at the bottom of the vial. Use this intermediate dilution when calculating final dilutions as recommended below. Store the vial at -20°C or below after dilution. Avoid cycles of freezing and thawing.

Note

For research use only

Application notes

Anti-Lysozyme Hen Egg White purified antibody has been tested for use in ELISA and western blot. Specific conditions for reactivity should be optimized by the end user. Expect a band at ~14.4-16.2 kDa in size corresponding to lysozyme by western blotting in the appropriate cell lysate or extract.

Isotype


IgG

Clonality


Polyclonal

Purity


This purified antibody is directed against lysozyme from hen egg white protein. This product is an IgG fraction antibody purified from monospecific antiserum by a multi-step process, which includes delipidation, salt fractionation and ion exchange chromatography followed by extensive dialysis against the buffer stated above. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Rabbit Serum as well as purified and partially purified Lysozyme (Hen Egg White). A BLAST analysis was used to suggest that this antibody would react with all forms of lysozyme from chicken, including



ELISA Results of Rabbit Anti-Lysozyme An...



Western Blot of Rabbit anti-Lysozyme (He...



Western Blot of Rabbit anti-Lysozyme ant...