

Product Datasheet

p21Cip1 (phospho-Ser130) antibody (orb34481)

Description	Rabbit polyclonal to p21Cip1 (phospho-Ser130)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT, WB
Immunogen	Synthetic Peptide
Target	CDKN1A
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	18119
Uniprot ID	P38936
NCBI	NP_000380.1 , NP_510867.1 , NP_001207706.1 , NP_001207707.1
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.

