




Product Datasheet

P21CIP1 (phospho-Thr145) antibody (orb34483)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT, IHC-P, WB
Immunogen	Synthetic Peptide
Target	CDKN1A
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	18119
Uniprot ID	P38936
NCBI	NP_000380.1 , NP_510867.1 , NP_001207706.1 , NP_001207707.1
Dilution Range	WB: 1:1000, DOT: 1:500, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Dot blot analysis of nitrocellulose memb...



Immunohistochemical analysis of formalin...