

Product Datasheet

p27Kip1 (phospho-Thr187) antibody (orb34488)

Description Peptide affinity purified rabbit polyclonal antibody.

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications IHC-P, WB

Immunogen Synthetic Peptide

Target CDKN1B

Form/Appearance Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal



MW 22073

Uniprot ID [P46527](#)

NCBI [NP_004055.1](#)

Dilution Range WB: 1:1000, IHC-P: 1:50-100

Expiration Date 12 months from date of receipt.



Western blot analysis of in HL60 tissue ...