

Product Datasheet

p53 (phospho-Ser315) antibody (orb34491)

Description

Rabbit polyclonal antibody to p53 (phospho-

P-Pab

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

DOT, IHC-P

Applications
Immunogen



Synthetic Peptide

Target

TP53

Form/Appearance

Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.



Immunohistochemical staining of paraffin...

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles

Note

For research use only

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

43653

Uniprot ID
P04637
NCBI
NP_001119588.1, NP_001263624.1, NP_001263627.1, NP_001119587.1, NP_001119589.1, NP_001263626.1, NP_001119586.1, NP_001119585.1, NP_000537.3, NP_00126, NP_001119590.1, NP_001263625.1, NP_001119584.1, NP_001263628.1
Dilution Range

DOT: 1:500, IHC-P: 1:50-100

Expiration Date

12 months from date of receipt.