

Product Datasheet

RAD9 (phospho-Ser277) antibody (orb34503)

Description Peptide affinity purified rabbit polyclonal antibody.

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications IHC-P, WB

Immunogen Synthetic Peptide

Target RAD9A

Form/Appearance Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal



MW 42547

Uniprot ID [Q99638](#)

NCBI [NP_001230153.1](#), [NP_004575.1](#)

Dilution Range WB: 1:1000, IHC-P: 1:50-100

Expiration Date 12 months from date of receipt.



Western blot analysis of in Y79 tissue l...