

Product Datasheet

MIP-3 alpha antibody (Peroxidase) (orb345230)

Description

MIP-3 alpha antibody (Peroxidase)

Species/Host

Rabbit

Reactivity

Mouse, Rat

Conjugation

HRP

Tested Applications

ELISA, IHC, WB

Immunogen

MIP-3 alpha Antibody Peroxidase
 Conjugated was prepared from whole rabbit serum produced by repeated immunizations with full length recombinant mouse MIP-3 α protein.

Preservatives

0.01% (w/v) Gentamicin Sulfate. Do NOT add Sodium Azide!

Form/Appearance

Lyophilized

Concentration

1.0 mg/mL

Storage

Store MIP-3 alpha Antibody Peroxidase Conjugated at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.

Note

For research use only

Application notes

This purified antibody has been tested in western blotting and is suitable for ELISA. By western blot a band approximately 10.7 kDa in size corresponding to mouse MIP-3 α protein is expected in the appropriate cell lysate or extract. Specific conditions for reactivity should be optimized by the end user.

Isotype

IgG

Clonality

Polyclonal

Purity

Mouse MIP-3 alpha Antibody Peroxidase Conjugated was heated to 56°C for 30 minutes. In ELISA and other immunoreactive assays, this antibody

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558