

---

## Product Datasheet

### GFP antibody (Biotin) (orb345335)

**Description**

GFP antibody (Biotin)

**Species/Host**

Mouse

**Reactivity**

Other

**Conjugation**

Biotin

**Tested**

ELISA, IHC, WB

**Applications**
**Immunogen**

Anti-Green Fluorescent Protein (GFP) is produced by immunizing GFP containing fusion protein that corresponds to the full length amino acid sequence (246aa) derived from the jellyfish *Aequorea victoria*.

**Preservatives**

0.01% (w/v) Sodium Azide

**Form/Appearance**

Lyophilized

**Concentration**

1.0 mg/mL

**Storage**

Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.

**Note**

For research use only

**Application notes**


Monoclonal anti-GFP is designed to detect enhanced GFP and GFP containing recombinant proteins. Tested in E, WB, IHC. This antibody can be used to detect GFP by ELISA (sandwich or capture) for the direct binding of antigen. Biotin conjugated monoclonal anti-GFP is well suited to titrate GFP in a sandwich ELISA in combination with Rockland's polyclonal anti-GFP (600-101-215) as the capture antibody. Only use the monoclonal form for the detection of enhanced or recombinant GFP. Polyclonal anti-GFP detects all variants of GFP tested to date. The biotin conjugated detection antibody is typically used with streptavidin conjugated HRP (code # S000-03) or other streptavidin conjugates. The use of polyclonal anti-GFP results in significant amplification of signal when fluorochrome conjugated polyclonal anti-GFP is used relative to the fluorescence of GFP alone. For immunoblotting use either alkaline phosphatase or peroxidase conjugated anti-GFP to detect GFP or GFP containing proteins on western blots. Optimal titers for applications should be determined by the researcher.

**Isotype**

IgG

**Clonality**

Monoclonal



IHC -  
Immuno-  
Fluorescence  
of Biotin  
Mous...



Western Blot  
of Anti-GFP  
(MOUSE)  
Monoclo...