



Product Datasheet

TrkA (phospho-Tyr791) antibody (orb34538)

Description	Rabbit polyclonal to TrkA (phospho-Tyr791)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT, WB
Immunogen	Synthetic Peptide
Target	NTRK1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	87497
Uniprot ID	P04629
NCBI	NP_002520.2 , NP_001007793.1 , NP_001012331.1
Dilution Range	WB: 1:1000, DOT: 1:500
Expiration Date	12 months from date of receipt.



WB analysis of A431 cell line using TrkA...