

## Product Datasheet

### FGFR1 (phospho-Tyr307) antibody (orb34568)

<b>Description</b>	Rabbit polyclonal to FGFR1 (phospho-Tyr307)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	DOT
<b>Immunogen</b>	Synthetic Peptide
<b>Target</b>	FGFR1
<b>Form/Appearance</b>	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
<b>Note</b>	For research use only
<b>Isotype</b>	Rabbit IgG
<b>Clonality</b>	Polyclonal
<b>MW</b>	91868
<b>Uniprot ID</b>	<a href="#">P11362</a>
<b>NCBI</b>	<a href="#">NP_001167538.1</a> , <a href="#">NP_001167536.1</a> , <a href="#">NP_075593.1</a> , <a href="#">NP_056934.2</a> , <a href="#">NP_075598.2</a> , <a href="#">NP_001167535.1</a> , <a href="#">NP_001167534.1</a> , <a href="#">NP_075594.1</a> , <a href="#">NP_001167537.1</a>
<b>Dilution Range</b>	DOT: 1:500
<b>Expiration Date</b>	12 months from date of receipt.

