

Product Datasheet

MAP3K1 (phospho-Thr1383) antibody (orb34578)

Description	Rabbit polyclonal to MAP3K1 (phospho-Thr1383)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT
Immunogen	Synthetic Peptide
Target	MAP3K1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	164470
Uniprot ID	Q13233
NCBI	NP_005912.1
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.

