

Product Datasheet

MEF2C (phospho-Thr20) antibody (orb34581)

Description	Rabbit polyclonal antibody to MEF2C (phospho-Thr20) belongs to the MEF2 family that is expressed in ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT
Immunogen	Synthetic Peptide
Target	MEF2C
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	51221
Uniprot ID	Q06413
NCBI	NP_001180277.1 , NP_002388.2 , NP_001180278.1 , NP_001180276.1 , NP_001180279.1 , NP_001124477.1
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.

