

Product Datasheet

GFP Western Blot Kit antibody (orb345938)

Description	GFP Western Blot Kit antibody
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	WB
Preservatives	Wash buffers MUST NOT contain SODIUM AZIDE or other inhibitors of peroxidase activity!
Form/Appearance	n/a
Storage	See GFP Western Blot kit insert for complete instructions.
Note	For research use only
Application notes	GFP Western Blot Kit is ideal for detection of GFP-tagged recombinant proteins by western blot. The GFP Western Blot kit makes use of Rockland's optimized anti-GFP antibody. The Chemiluminescent GFP Western Blot kit is useful for both western blotting and dot blotting methods. Please read the entire product insert prior to use.
Clonality	Monoclonal
Purity	Biorbyt' Chemiluminescent Western Blot Kit for GFP combines all of the necessary reagents with a rapid proven protocol and the extremely high signal detection of our luminol chemiluminescent substrate for the detection of recombinant proteins containing GFP and its variants. The Chemiluminescent Western Blot Kit design includes straightforward procedures and color-coding to allow for ease of use. This kit contains: GFP protein, Mouse IgG control, mouse anti-GFP antibody, anti-mouse IgG peroxidase, Tween-20, BSA, and Femtomax, along with additional supplies and protocol; sufficient substrate for up to 30 mini blots at 7.5 x 8 cm ² (1,800 cm ²) and is stable for at least 1 year when stored as indicated.
Dilution Range	WB: User Optimized
Expiration Date	Please enquire.



Flow Diagram for Chemiluminescent Western...



Western Blot of anti-GFP monoclonal anti...