

Product Datasheet

Bovine IgG F(ab')₂ Agarose Antibody (orb345988)

Description	Bovine IgG F(ab') ₂ Agarose Antibody
Conjugation	Agarose
Tested Applications	IP
Preservatives	0.01% (w/v) Sodium Azide
Form/Appearance	Suspension of agarose beads
Concentration	0.5 mg/cc agarose
Storage	Store vial at 4° C prior to opening. DO NOT FREEZE.
Note	For research use only
Application notes	Bovine IgG F(ab') ₂ Agarose is suitable for immunoprecipitation and Immunoaffinity purifications.
Isotype	IgG F(ab') ₂
Purity	This product is purified Bovine IgG F(ab') ₂ isolated from normal Bovine Serum and then digested by the enzyme pepsin followed by subsequent chromatographic separation and coupling to activated agarose. Sufficient antibody capacity is provided to bind a minimum of 5.0 mg of anti-Bovine IgG.
Source	Bovine
Dilution Range	IP: User Optimized
Expiration Date	12 months from date of receipt.