

Product Datasheet

RPS6KB1 (phospho-Ser418) antibody (orb34607)

Description	Rabbit polyclonal to RPS6KB1 (phospho-Ser418)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT
Immunogen	Synthetic Peptide
Target	RPS6KB1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	59140
Uniprot ID	P23443
NCBI	NP_003152.1 , NP_001258989.1
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.

