

Product Datasheet

Goat IgG Agarose Antibody (orb346104)

Description	Goat IgG Agarose Antibody
Conjugation	Agarose
Tested Applications	SDS-PAGE
Preservatives	0.05% (w/v) Sodium Azide
Form/Appearance	Suspension of agarose beads
Concentration	8.56 mg/cc packed agarose
Storage	Store vial at 4° C prior to opening. DO NOT FREEZE.
Note	For research use only
Application notes	Goat IgG Agarose has been tested in SDS-Page and is suitable for immunoprecipitation and Immunoaffinity purifications.
Isotype	IgG
Purity	This product is normal Goat IgG coupled to activated agarose. A single precipitin arc was observed against anti-Goat Serum when assayed by immunoelectrophoresis prior to coupling to the beads.
Source	Goat
Expiration Date	12 months from date of receipt.



SDS-PAGE of
Goat IgG
Whole
Molecule
Rhod...