


## Product Datasheet

### Goat IgG Fab Peroxidase Antibody (orb346118)

<b>Description</b>	Goat IgG Fab Peroxidase Antibody
<b>Conjugation</b>	HRP
<b>Tested Applications</b>	SDS-PAGE
<b>Preservatives</b>	0.01% (w/v) Gentamicin Sulfate. Do NOT add Sodium Azide!
<b>Form/Appearance</b>	Lyophilized
<b>Concentration</b>	1.0 mg/mL
<b>Storage</b>	Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
<b>Note</b>	For research use only
<b>Application notes</b>	Goat IgG F(ab) fragment peroxidase conjugate has been tested by SDS-PAGE and can be utilized as a control or standard reagent in Western Blotting and ELISA experiments.
<b>Isotype</b>	IgG Fab
<b>Purity</b>	This product was prepared from normal serum by delipidation, salt fractionation, ion exchange chromatography followed by papain digestion and extensive dialysis against the buffer stated above. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Peroxidase, anti-Goat IgG, anti-Goat IgG F(ab') <sub>2</sub> and anti-Goat Serum. No reaction was observed against anti-Goat IgG F(c) or anti-Papain.
<b>Source</b>	Goat
<b>Expiration Date</b>	12 months from date of receipt.



SDS-PAGE of  
Goat IgG  
Whole  
Molecule  
Rhod...