

Product Datasheet

ATF2 (phospho-Ser322) antibody (orb34619)

Description	Rabbit polyclonal to ATF2 (phospho-Ser322)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT
Immunogen	Synthetic Peptide
Target	ATF2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	54537
Uniprot ID	P15336
NCBI	NP_001243020.1 , NP_001871.2 , NP_001243019.1 , NP_001243023.1 , NP_001243022.1 , NP_001243021.1
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.

