

## Product Datasheet

### Human IgG Agarose Antibody (orb346215)

<b>Description</b>	Human IgG Agarose Antibody
<b>Conjugation</b>	Agarose
<b>Tested Applications</b>	SDS-PAGE
<b>Preservatives</b>	0.01% (w/v) Sodium Azide
<b>Form/Appearance</b>	Suspension of agarose beads
<b>Concentration</b>	2.5 mg/mL
<b>Storage</b>	Store vial at 4° C prior to opening. DO NOT FREEZE.
<b>Note</b>	For research use only
<b>Application notes</b>	Human IgG Agarose has been tested by SDS-Page and is suitable for immunoprecipitation and Immunoaffinity purifications.
<b>Isotype</b>	IgG
<b>Purity</b>	This product was prepared by coupling purified Human IgG antibody to activated agarose. Sufficient antibody capacity is provided to bind a minimum of 5 mg of anti-Human IgG antibody.
<b>Source</b>	Human
<b>Expiration Date</b>	12 months from date of receipt.