


Product Datasheet

Rabbit IgG Agarose Antibody (orb346307)

Description	Rabbit IgG Agarose Antibody
Conjugation	Agarose
Tested Applications	IP
Preservatives	0.01% (w/v) Sodium Azide
Form/Appearance	Suspension of agarose beads
Concentration	10 mg/mL
Storage	Store vial at 4° C prior to opening. DO NOT FREEZE.
Note	For research use only
Application notes	Rabbit IgG Agarose is suitable for immunoprecipitation and Immunoaffinity purifications.
Isotype	IgG
Purity	This product is normal Rabbit IgG coupled to activated agarose. A single precipitin arc was observed against anti-Rabbit Serum when assayed by immunoelectrophoresis prior to coupling to the beads.
Source	Rabbit
Dilution Range	IP: User Optimized
Expiration Date	12 months from date of receipt.



SDS-PAGE of Rabbit IgG Whole Molecule Rh...