


Product Datasheet

Rabbit IgG F(ab')₂ Agarose Antibody (orb346308)

Description	Rabbit IgG F(ab') ₂ Agarose Antibody
Conjugation	Agarose
Tested Applications	IP
Preservatives	0.5% (w/v) Sodium Azide
Form/Appearance	Suspension of agarose beads
Storage	Store vial at 4° C prior to opening. DO NOT FREEZE.
Note	For research use only
Application notes	Rabbit IgG F(ab') ₂ Agarose is suitable for immunoprecipitation and Immunoaffinity purifications.
Isotype	IgG F(ab') ₂
Purity	This product is an purified Rabbit IgG isolated from normal Rabbit Serum digested by the enzyme pepsin and purified by chromatography coupled to activated agarose. Sufficient antibody capacity is provided to bind a minimum of 5 mg of anti-Rabbit IgG.
Source	Rabbit
Dilution Range	IP: User Optimized
Expiration Date	12 months from date of receipt.



SDS-PAGE of Rabbit IgG Whole Molecule Rh...