

Product Datasheet

MEF2C (phospho-Ser387) antibody (orb34631)

Description Rabbit polyclonal antibody to MEF2C (phospho-Thr20) belongs to the MEF2 family that is expressed in ...

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications DOT

Immunogen Synthetic Peptide

Target MEF2C

Form/Appearance Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal

MW 51221

Uniprot ID [Q06413](#)

NCBI [NP_001180277.1](#), [NP_002388.2](#), [NP_001180278.1](#), [NP_001180276.1](#), [NP_001180279.1](#), [NP_001124477.1](#)

Dilution Range DOT: 1:500

Expiration Date 12 months from date of receipt.

