


Product Datasheet

Rat IgG Agarose Antibody (orb346335)

Description	Rat IgG Agarose Antibody
Conjugation	Agarose
Tested Applications	IP
Preservatives	0.01% (w/v) Sodium Azide
Form/Appearance	Suspension of agarose beads
Storage	Store vial at 4° C prior to opening. DO NOT FREEZE.
Note	For research use only
Application notes	Rat IgG Agarose is suitable for immunoprecipitation and Immunoaffinity purifications.
Isotype	IgG
Purity	This product is normal Rat IgG coupled to activated agarose. A single precipitin arc was observed against anti-Rat Serum when assayed by immunoelectrophoresis prior to coupling to the beads. Sufficient capacity is provided to bind a minimum of 10 mg of anti-Rat IgG.
Source	Rat
Dilution Range	IP: User Optimized
Expiration Date	12 months from date of receipt.



SDS-PAGE of Rat IgG F(c) Fragment Rhodam...