


Product Datasheet

Rat IgG F(ab')₂ Agarose Antibody (orb346336)

Description	Rat IgG F(ab') ₂ Agarose Antibody
Conjugation	Agarose
Tested Applications	IP
Preservatives	0.01% (w/v) Sodium Azide
Form/Appearance	Suspension of agarose beads
Storage	Store vial at 4° C prior to opening. DO NOT FREEZE.
Note	For research use only
Application notes	Rat IgG F(ab') ₂ Agarose is suitable for immunoprecipitation and Immunoaffinity purifications.
Isotype	IgG F(ab') ₂
Purity	This product is an purified Rat IgG isolated from normal Rat Serum digested by the enzyme pepsin and purified by chromatography coupled to activated agarose. Sufficient antibody capacity is provided to bind a minimum of 5.0 of anti-Rat IgG.
Source	Rat
Dilution Range	IP: User Optimized
Expiration Date	12 months from date of receipt.



SDS-PAGE of Rat IgG F(c) Fragment Rhodam...