

Product Datasheet

Bovine IgG antibody (Agarose) (orb346548)

Description	Bovine IgG antibody (Agarose)
Species/Host	Goat
Reactivity	Bovine
Conjugation	Agarose
Tested Applications	SDS-PAGE
Preservatives	0.01% (w/v) Sodium Azide
Form/Appearance	Suspension of agarose beads
Concentration	9.2 mg/cc
Storage	Store vial at 4° C prior to opening. DO NOT FREEZE.
Note	For research use only
Application notes	Anti-Bovine IgG [H&L] (GOAT) Antibody Agarose Conjugated has been tested by SDS-Page and is applicable to purification of Bovine IgG from sample solution.
Isotype	IgG
Clonality	Polyclonal
Purity	Anti-Bovine IgG [H&L] (GOAT) Antibody Agarose Conjugated is an IgG fraction antibody coupled to activated agarose. Sufficient antibody capacity is provided to bind a minimum of 5 mg of pure Bovine IgG.
Expiration Date	12 months from date of receipt.