

Product Datasheet

Goat IgG antibody (Agarose) (orb346594)

Description	Goat IgG antibody (Agarose)
Species/Host	Rabbit
Reactivity	Goat
Conjugation	Agarose
Preservatives	0.01% (w/v) Sodium Azide
Form/Appearance	Suspension of agarose beads
Concentration	10 mg/mL
Storage	Store Anti-Goat IgG [H&L] (Goat) Antibody Agarose Conjugated vial at 4° C prior to opening. DO NOT FREEZE.
Note	For research use only
Application notes	Anti-Goat IgG [H&L] (Goat) Antibody Agarose Conjugated is applicable to purification of Goat IgG from sample solution.
Isotype	IgG
Clonality	Polyclonal
Purity	Anti-Goat IgG [H&L] (Rabbit) Antibody Agarose Conjugated is an IgG fraction antibody coupled to activated agarose. Sufficient antibody capacity is provided to bind a minimum of 5 mg of pure Goat IgG.
Expiration Date	12 months from date of receipt.