
Product Datasheet

Human IgG F(ab')₂ antibody (FITC) (orb346658)

Description

 Human IgG F(ab')₂ antibody (FITC)

Species/Host

Rabbit

Reactivity

Human

Conjugation

FITC

Tested Applications

FC, FLISA, IF

Immunogen

 Human IgG F(ab')₂ fragment

Preservatives

0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2

Form/Appearance

Lyophilized

Concentration

10 mg/mL

Storage

Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.

Note

For research use only

Application notes

Anti-Human IgG F(ab')₂ Fluorescein Antibody is designed for immunofluorescence microscopy, fluorescence based plate assays (FLISA) and fluorescent western blotting. This product is also suitable for multiplex analysis, including multicolor imaging, utilizing various commercial platforms.

Isotype

IgG

Clonality

Polyclonal

Purity

This product is an IgG fraction antibody purified from monospecific antiserum by a multi-step process which includes delipidation, salt fractionation and ion exchange chromatography followed by extensive dialysis against the buffer stated above. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-fluorescein, anti-Rabbit Serum, Human IgG, Human IgG F(ab')₂ and Human Serum. No reaction was observed against