
Product Datasheet

Mouse IgM (H&L) antibody (orb346666)

Analysis Reagents
Description

Mouse IgM (H&L) antibody

Species/Host

Goat

Reactivity

Mouse

Conjugation

Unconjugated

Tested Applications

ELISA, IHC, WB

Immunogen

Mouse IgM whole molecule

Preservatives

0.01% (w/v) Sodium Azide

Form/Appearance

Lyophilized

Concentration

10.0 mg/mL

Storage

Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.

Note

For research use only

Application notes

Anti-Mouse IgM (mu chain) (GOAT) Antibody has been tested by ELISA and is assayed against Mouse IgM in a standard capture ELISA using Peroxidase Conjugated Streptavidin. A working dilution of 1:15,000 to 1:70,000 of the concentration is suggested for this product. Mouse IgM antibody is suitable for immunoassays including western blot and ELISA.

Isotype

IgG

Clonality

Polyclonal

Purity

Anti-MOUSE IgM (H&L) (GOAT) Antibody is an IgG fraction antibody purified from monospecific antiserum by a multi-step process which includes delipidation, salt fractionation and ion exchange chromatography followed by extensive dialysis against the buffer stated above. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Goat Serum, Mouse IgM and Mouse Serum. Some light chain cross reactivity may occur against Mouse IgG.