

Product Datasheet

RABBIT IgG antibody (Agarose) (orb346685)

Description	RABBIT IgG antibody (Agarose)
Species/Host	Goat
Reactivity	Rabbit
Conjugation	Agarose
Tested Applications	IP
Preservatives	0.1% (w/v) Sodium Azide
Form/Appearance	Suspension of agarose beads
Concentration	0.5 cc agarose per ml of suspension
Storage	Store vial at 4° C prior to opening. DO NOT FREEZE.
Note	For research use only
Application notes	Anti Rabbit IgG Secondary Antibody Agarose can be used for immunoaffinity chromatography of polyclonal or monoclonal rabbit IgG antibodies. Once bound, non-specific antibodies are removed by copious washing followed by elution of specific antibody at low pH.
Isotype	IgG
Clonality	Polyclonal
Purity	This product is an IgG fraction antibody coupled to activated agarose. Sufficient antibody capacity is provided to bind a minimum of 5 mg of pure Rabbit IgG.
Dilution Range	IP: User Optimized
Expiration Date	12 months from date of receipt.