
Product Datasheet

**F(ab')₂ Chicken IgG (H&L) antibody (Peroxidase)
(orb346746)**

Description

 F(ab')₂ Chicken IgG (H&L) antibody

Species/Host	Rabbit
Reactivity	Gallus
Conjugation	HRP
Tested Applications	ELISA, IHC, WB
Immunogen	Chicken IgG whole molecule
Preservatives	0.01% (w/v) Gentamicin Sulfate. Do NOT add Sodium Azide!
Form/Appearance	Lyophilized
Concentration	10.0 mg/mL
Storage	Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Note	For research use only
Application notes	F(ab') ₂ Chicken IgG Antibody Peroxidase Conjugated is ideal for investigators who routinely perform flow cytometry, immunohistochemistry or IHC and other immunoassays.
Isotype	IgG F(ab') ₂
Clonality	Polyclonal
Purity	This product is a F(ab') ₂ fragment of an IgG fraction antibody purified from monospecific antiserum by a multi-step process which includes delipidation, salt fractionation, ion exchange chromatography and pepsin digestion followed by extensive dialysis against the buffer stated above. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-