
Product Datasheet

**F(ab')₂ Golden Syrian Hamster IgG (H&L) antibody
(FITC) (orb346765)**

Description

 F(ab')₂ Golden Syrian Hamster IgG

Species/Host	Rabbit
Reactivity	Golden Hamster
Conjugation	FITC
Tested Applications	FC, FLISA, IF
Immunogen	Hamster IgG whole molecule
Preservatives	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Form/Appearance	Lyophilized
Storage	Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Note	For research use only
Application notes	This product is designed for immunofluorescence microscopy, fluorescence based plate assays (FLISA) and fluorescent western blotting. This product is also suitable for multiplex analysis, including multicolor imaging, utilizing various commercial platforms.
Isotype	IgG F(ab') ₂
Clonality	Polyclonal
Purity	This product is a F(ab') ₂ fragment of an IgG fraction antibody purified from monospecific antiserum by a multi-step process which includes delipidation, salt fractionation, ion exchange chromatography and pepsin digestion followed by extensive dialysis against the buffer stated above. Assay by immunoelectrophoresis resulted in a

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

 Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709. United States

 Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Dilution Range

FLISA: 1:10,000 - 1:50,000. FC: