
Product Datasheet

Chicken IgG (H&L) antibody (Biotin) (orb346896)

Description

Chicken IgG (H&L) antibody (Biotin)

Species/Host

Goat

Reactivity


Gallus

Conjugation

Biotin

Tested

ELISA, IHC, WB

Applications

 ELISA
results of
purified
Goat anti-
Chic...

Immunogen

Chicken IgG whole molecule

Preservatives

0.01% (w/v) Sodium Azide

Form/Appearance

Lyophilized

Concentration

2.0 mg/mL

Storage

Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.

Note

For research use only

Application notes

Anti-Chicken IgG Biotin conjugated antibody has been tested by ELISA and is suitable for ELISA, Western Blot and Immunohistochemistry applications. Specific conditions for reactivity should be optimized by end user.

Isotype

IgG

Clonality

Polyclonal

Purity

Anti-Chicken IgG antibody was prepared from monospecific antiserum by immunoaffinity chromatography using Chicken IgG coupled to agarose beads. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-biotin, anti-Goat Serum, Chicken IgG and Chicken Serum.

Dilution Range

ELISA: 1:98,500, IHC: 1:500 - 1:3,000, WB: 1:3,000 - 1:30,000

Expiration Date

12 months from date of receipt.