

Product Datasheet

PDPK1 (phospho-Ser396) antibody (orb34701)

Description	Rabbit polyclonal to PDPK1 (phospho-Ser396)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT
Immunogen	Synthetic Peptide
Target	PDPK1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	63152
Uniprot ID	O15530
NCBI	NP_002604.1 , NP_112558.2 , NP_001248745.1
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.

