

---

## Product Datasheet

### Goat IgG (H&L) antibody (Texas Red) (orb347020)

Analysis Reagents

Goat IgG (H&L) antibody (Texas Red)

<b>Description</b>	Goat IgG (H&L) antibody (Texas Red)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Goat
<b>Conjugation</b>	Texas Red
<b>Tested Applications</b>	FLISA, IF, WB
<b>Immunogen</b>	Goat IgG whole molecule
<b>Preservatives</b>	0.01% (w/v) Sodium Azide
<b>Form/Appearance</b>	Lyophilized
<b>Concentration</b>	2.0 mg/mL
<b>Storage</b>	Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
<b>Note</b>	For research use only
<b>Application notes</b>	Anti-Goat IgG Texas Red Antibody has been tested by Western blot. This product is designed for immunofluorescence microscopy, fluorescence based plate assays (FLISA) and fluorescent western blotting. This product is also suitable for multiplex analysis, including multicolor imaging, utilizing various commercial platforms.
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Purity</b>	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Goat IgG coupled to agarose. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Rabbit Serum, Goat IgG and Goat Serum.
<b>Dilution Range</b>	FLISA: 1:10,000 - 1:50,000, IF: 1:1,000 - 1:5,000
<b>Expiration Date</b>	12 months from date of receipt.