

Product Datasheet

PIK3CG (phospho-Ser1100) antibody (orb34707)

Description	Rabbit polyclonal to PIK3CG (phospho-Ser1100)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT
Immunogen	Synthetic Peptide
Target	PIK3CG
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	126454
Uniprot ID	P48736
NCBI	NP_002640.2 , NP_001269355.1 , NP_001269356.1
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.

