



Product Datasheet

ABL (phospho-Tyr204) antibody (orb34723)

Description	Rabbit polyclonal to ABL (phospho-Tyr204)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT, IF
Immunogen	Synthetic Peptide
Target	ABL1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	122873
Uniprot ID	P00519
NCBI	NP_005148.2 , NP_009297.2
Dilution Range	DOT: 1:500, IF/ICC: 1:10-50
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of MCF-7 cel...