

Product Datasheet

MBP (phospho-Tyr203) antibody (orb34793)

Description	Rabbit polyclonal to MBP (phospho-Tyr203)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT
Immunogen	Synthetic Peptide
Target	MBP
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	33117
Uniprot ID	P02686
NCBI	NP_002376.1 , NP_001020261.1 , NP_001020252.1 , NP_001020271.1 , NP_001020272.1 , NP_001020263.1
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.

