
Product Datasheet

Sheep IgG (H&L) antibody (Texas Red) (orb347943)

Evolution Bioproducts.

Sheep IgG (H&L) antibody (Texas Red)

Description	Sheep IgG (H&L) antibody (Texas Red)
Species/Host	Rabbit
Reactivity	Sheep
Conjugation	Texas Red
Tested Applications	DOT, FC, FLISA, IF, WB
Immunogen	Sheep IgG whole molecule
Preservatives	0.01% (w/v) Sodium Azide
Form/Appearance	Lyophilized
Concentration	1.0 mg/mL
Storage	Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Note	For research use only
Application notes	Anti-Sheep IgG (H&L) Texas Red Antibody has been tested by dot blot and western blot and is designed for immunofluorescence microscopy, fluorescence based plate assays (FLISA) and fluorescent western blotting. This product is also suitable for multiplex analysis, including multicolor imaging, utilizing various commercial platforms.
Isotype	IgG
Clonality	Polyclonal
Purity	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Sheep IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Rabbit Serum, Sheep IgG and Sheep Serum.
Dilution Range	FLISA: 1:10,000 - 1:50,000, FC: 1:500 -