

Product Datasheet

GFAP (phospho-Ser8) antibody (orb34799)

Description Rabbit polyclonal antibody to GFAP is a member of the intermediate filament family that is expressed...

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications DOT

Immunogen Synthetic Peptide

Target GFAP

Form/Appearance Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal

MW 49880

Uniprot ID [P14136](#)

NCBI [NP_002046.1](#), [NP_001229305.1](#), [NP_001124491.1](#)

Dilution Range DOT: 1:500

Expiration Date 12 months from date of receipt.

P-Pab

