

---

## Product Datasheet

**Monkey IgG IgA IgM (H&L) antibody (FITC)  
(orb348016)**

**Description**

Monkey IgG IgA IgM (H&amp;L) antibody (FITC)

**Species/Host**

Goat

**Reactivity**

Monkey

**Conjugation**

FITC

**Tested**

FC, FLISA, IF, WB

**Applications**
**Immunogen**

Monkey IgG, IgA and IgM whole molecules

**Preservatives**

0.01% (w/v) Sodium Azide

**Form/Appearance**

Lyophilized

**Concentration**

1.0 mg/mL

**Storage**

Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.

**Note**

For research use only

**Application notes**

Anti-Monkey IgG IgA IgM Fluorescein Conjugated Antibody has been tested by western blot and is designed for immunofluorescence microscopy, fluorescence based plate assays (FLISA) and fluorescent western blotting. This product is also suitable for multiplex analysis, including multicolor imaging, utilizing various commercial platforms.

**Isotype**

IgG

**Clonality**

Polyclonal

**Purity**

This product was prepared from polyspecific antiserum by immunoaffinity chromatography using antigens coupled to agarose beads. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Fluorescein, anti-Goat Serum, Monkey IgG, Monkey IgA and Monkey IgM. This reagent is suitable for the detection of all Monkey isotypes and chain combinations.

**Dilution Range**

FLISA: 1:10,000 - 1:50,000, FC: 1:500 - 1:2,500, IF: 1:1,000 - 1:5,000

**Expiration Date**

12 months from date of receipt.

Product	Reactivity	Reactivity	Reactivity
Fluorescein	+	+	+
Fluorescein	+	+	+
Fluorescein	+	+	+
Fluorescein	+	+	+
Fluorescein	+	+	+
Fluorescein	+	+	+

This table displays additional reactivit...