
Product Datasheet

Golden Syrian & Armenian Hamster IgG (H&L) antibody (Peroxidase) (orb348069)

Description

Golden Syrian & Armenian

(Peroxidase)

Species/Host

Goat

Reactivity

Golden Hamster, Hamster

Conjugation

HRP

Tested Applications

ELISA, IHC, WB

Immunogen

Armenian and Golden Syrian Hamster IgG, whole molecule

Preservatives

 0.01% (w/v) Gentamicin Sulfate.
Do NOT add Sodium Azide!

Form/Appearance

Lyophilized

Concentration

1.0 mg/mL

Storage

Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.

Note

For research use only

Application notes

Anti-Golden Syrian & Armenian Hamster IgG Peroxidase conjugate is suitable for immunoblotting (western or dot blot), ELISA, immunoelectron microscopy and immunohistochemistry as well as other antibody-based enzymatic assays requiring lot-to-lot consistency.

Isotype

IgG

Clonality

Polyclonal

Purity

This product was prepared from monospecific antiserum by immunoaffinity chromatography using Armenian/Golden Syrian Hamster IgG coupled to agarose