

Product Datasheet

MDM2 (phospho-Ser395) antibody (orb34816)

Description	Rabbit polyclonal antibody to MDM2 (phospho-Ser395) also known as HDM2, MTBP, p53 Binding Protein, b...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT
Immunogen	Synthetic Peptide
Target	MDM2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	55233
Uniprot ID	Q00987
NCBI	NP_001265391.1 , NP_001138811.1 , NP_002383.2
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.

