

Product Datasheet

PLXND1 (phospho-Tyr1642) antibody (orb34821)

Description	Rabbit polyclonal to PLXND1 (phospho-Tyr1642)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT
Immunogen	Synthetic Peptide
Target	PLXND1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	212007
Uniprot ID	Q9Y4D7
NCBI	NP_055918.2
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.

