

Product Datasheet

ANTXR1 (phospho-Tyr425) antibody (orb34840)

Description	Rabbit polyclonal to ANTXR1 (phospho-Tyr425)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT
Immunogen	Synthetic Peptide
Target	ANTXR1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	62789
Uniprot ID	Q9H6X2
NCBI	NP_060623.2 , NP_115584.1 , NP_444262.1
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.

