


Product Datasheet

Rat PD-1 / PDCD1 Protein (orb348794)

Description	Rat PD-1, His Tag is expressed from human 293 cells (HEK293). It contains AA Leu 25 - Gln 167 (Acces...
Reactivity	Rat
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	PD-1
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 18.1 kDa. The protein migrates as 30-44 kDa under reducing (R) condition (SDS-PAGE).
Protein Sequence	NP_001100397.1
Purity	95%
Source	Rat PD-1, His Tag (orb348794) is expressed from human 293 cells (HEK293). It contains AA Leu 25 - Gln 167 (Accession # NP_001100397.1).
MW	18.1 kDa
Uniprot ID	D3ZIN8
NCBI	NP_001100397.1
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Rat PD-1 protein