


Product Datasheet

Mouse BTLA Protein (orb348804)

Description	Mouse BTLA, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Glu 30 - Gly 176 (Acce...
Reactivity	Mouse
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	BTLA
Preservatives	Tris with Glycine, Arginine and NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 42.8 kDa. The protein migrates as 55-66 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	Glu 30 - Gly 176
Purity	95%
Source	Mouse BTLA, Fc Tag (orb348804) is expressed from human 293 cells (HEK293). It contains AA Glu 30 - Gly 176 (Accession # Q32MV9).
MW	42.8 kDa
Uniprot ID	Q32MV9
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Mouse BTLA protein