


Product Datasheet

Human FKBP4 / FKBP52 Protein (orb348809)

Description	Human FKBP4, His,PA Tag is expressed from human 293 cells (HEK293). It contains AA Met 1 - Ala 459 (...)
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	FKBP4
Preservatives	50 mM Tris, 100 mM Glycine, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His & C-PA Tag
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus, followed by a PA tag, and has a calculated MW of 55.1 kDa. The reducing (R) protein migrates as 60 kDa in SDS-PAGE.
Protein Sequence	AAH01786
Purity	95%
Source	Human FKBP4, His,PA Tag (orb348809) is expressed from human 293 cells (HEK293). It contains AA Met 1 - Ala 459 (Accession # AAH01786).
MW	55.1 kDa
Uniprot ID	Q02790
NCBI	AAH01786
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human FKBP4 protein